

Endomotor



Economical .

Torque Control.

Auto Reverse.



Customer

Service

Support

Email : healixhealthcare@gmail.com

Tel. : +919324607948 / +919223444316



Endomotor

Economical. Torque Control. Auto Reverse.



Speed Range*:
126 to 625 rpm

Torque Range*:
3 to 52mNm

Auto Reverse
Auto Stop
Auto 1 cycle reverse

9 adjustable
programs

Works with all
rotary file systems

How to Connect the Parts:

- ▶ Connect the Foot Switch to the Foot Switch Connector on the Control Panel
- ▶ Connect the Motor to the Motor Connector on the Control Panel
- ▶ Plug the Power Cable into the Power Socket of 230V

Parts

1. Control Panel
2. Motor
3. Foot Switch

110, Manish Ind. Estate. No. 4, Navghar Road, Vasai East), Dist: Palghar - 401210, India, Tel: 0250-2390707 / 2393973

www.dentaldeal.in | Email : healixhealthcare@gmail.com | Tel. : +91 9324 6079 48 / +91 9223 4443 16